

Kallikrein 14 Antibody

Catalog No: #39246

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Kallikrein 14 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified
Applications	ELISA WB
Species Reactivity	Hu
Specificity	Kallikrein 14 antibody detects endogenous levels of Kallikrein 14 protein. Not tested with other proteins
Target Name	Kallikrein 14
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA WB: 1:100

Images



Note: This product is for in vitro research use only