

## PP2A alpha and beta Antibody

Catalog No: #39248

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PP2A alpha and beta Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified
Applications	ELISA WB
Species Reactivity	Hu Ms Rt
Specificity	PP2A alpha and beta Antibody detects endogenous levels of PP2A alpha and beta protein. Not tested with other proteins
Immunogen Type	Peptide-KLH
Immunogen Description	Synthetic peptide conjugated to KLH
Target Name	PP2A alpha and beta
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

## Application Details

ELISA WB: 1:250

## Images



Note: This product is for in vitro research use only