

HES-1 Antibody

Catalog No: #39283

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HES-1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Not yet tested in other species.
Immunogen Type	Peptide-KLH
Immunogen Description	15 amino acid peptide from near the N terminal residues of human HES-1 protein. The peptide sequence is: CDALKKDSSRHSKLEK
Target Name	HES-1
Calculated MW	30kD
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

WB :1:1000-3000

Images



Western blot analysis Hes1 using 293 whole cell lysates

Note: This product is for in vitro research use only