

PLA2G12B Antibody

Catalog No: #39285

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PLA2G12B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	purified
Applications	ELISA WB
Species Reactivity	Hu Ms Rt
Specificity	Not yet tested in other species.
Immunogen Type	Peptide-KLH
Immunogen Description	a 15 amino acid peptide from near the N terminal residues of human PLA2G12B protein.
Target Name	PLA2G12B
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA WB

Images



Note: This product is for in vitro research use only