

NFkB3 Antibody

Catalog No: #39294

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NFkB3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Hu Rt Ms
Specificity	NFkB3 antibody detects endogenous levels of NFkB3 protein. Not tested with other proteins
Immunogen Type	Peptide-KLH
Immunogen Description	Synthetic peptide conjugated to KLH
Target Name	NFkB3
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

ELISA WB

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded amygdala showing nuclear staining with anti-NFkB3 antibody dilution: 1:100



Note: This product is for in vitro research use only