

## Stilbene Synthase Antibody

Catalog No: #39305

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Stilbene Synthase Antibody
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Vitis riparia
Immunogen Type	Peptide-KLH
Immunogen Description	peptide from near the C terminal residues of Vitis riparia stilbene synthase protein.
Target Name	Stilbene Synthase
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

## Application Details

ELISA: 1:10000-1:80000 WB: 1:200-1:1000

Note: This product is for in vitro research use only