

## AICDA Antibody

Catalog No: #39314

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	AICDA Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Ms
Immunogen Type	Peptide-KLH
Immunogen Description	Rabbit polyclonal AICDA(1) antibody was raised against a 16 amino acid peptide corresponding to near the N terminal residues of Human AICDA.
Target Name	AICDA
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


## Application Details

ELISA WB IHC

## Images



Immunohistochemical staining of formalin-fixed paraffin-embedded lymph node showing Cytoplasmic staining with anti-AICDA antibody dilution: 1:100.



---

Note: This product is for in vitro research use only