

Cyclin D1 Antibody

Catalog No: #39315

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Cyclin D1 Antibody
Host Species	rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide-KLH
Immunogen Description	Rabbit polyclonal Cyclin D1 antibody was raised against a 12 amino acid peptide from the C terminal residues of human Cyclin D1 protein.
Target Name	Cyclin D1
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

ELISA WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded brain showing nuclear staining with anti-Cyclin D1 antibody dilution: 1:100.



Note: This product is for in vitro research use only