

TPM4 Antibody

Catalog No: #39336

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TPM4 Antibody
Host Species	Rabbit
Clonality	PolyClonal
Purification	Affinity purified
Applications	ELISA IHC
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide-KLH
Immunogen Description	Rabbit polyclonal TPM4 antibody was raised against a 16 amino acid peptide from near the N terminal residues of human TPM4 protein.
Target Name	TPM4
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded skeletal muscle showing cytoplasmic staining with anti-TPM4 antibody dilution: 1:100.

Note: This product is for in vitro research use only