

TLE4 Antibody

Catalog No: #39338

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TLE4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA IHC WB
Species Reactivity	Hu Ms Rt
Specificity	The species of the brain tissue is Human. The peptide sequence of this antibody is corresponding to amino acids of human 1-50aa. This antibody is predicted to recognize mouse and rattus TLE4 protein because of synthesized peptide amino acid 100% homology with mouse and rattus TLE4 protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Rabbit polyclonal TLE4 antibody was raised against a peptide from N terminal residues of human TLE4 protein.
Target Name	TLE4
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

ELISA IHC WB

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded human brain showing nuclear staining with anti-TLE4 antibody dilution: 1:200.



Note: This product is for in vitro research use only