

ACOT7 Antibody

Catalog No: #39350

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ACOT7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA IHC WB
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide-KLH
Immunogen Description	Rabbit polyclonal ACOT7 antibody was raised against a peptide near the C terminal residues of Human ACOT7.
Target Name	ACOT7
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

ELISA IHC WB

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded brain showing cytoplasmic staining with anti-ACOT7 antibody at a dilution of 1/100.



Note: This product is for in vitro research use only