

ADORA1 Antibody

Catalog No: #39355

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ADORA1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA IHC
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide-KLH
Immunogen Description	Rabbit polyclonal ADORA1 antibody was raised against a peptide from intermediate residues of human ADORA1.
Target Name	ADORA1
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded rat kidney showing membrane staining with anti-ADORA1 antibody at a dilution of 1/100.

Note: This product is for in vitro research use only