

CIDE B Antibody

Catalog No: #39362

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | CIDE B Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Affinity purified |
| Applications | ELISA IHC |
| Species Reactivity | Hu |
| Immunogen Type | Peptide-KLH |
| Immunogen Description | Rabbit polyclonal CIDE B antibody was raised against a 14 amino acid peptide at the C terminal residues of human CIDE B protein. |
| Target Name | CIDE B |
| Concentration | 1.0mg/ml |
| Formulation | PBS, pH 7.4 with 0.02% Sodium Azide |
| Storage | This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt. |

Application Details

ELISA IHC

Note: This product is for in vitro research use only