

## N-acetylglucosaminyltransferase V Antibody

Catalog No: #39378

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	N-acetylglucosaminyltransferase V Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA IHC
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide-KLH
Immunogen Description	Rabbit polyclonal N-acetylglucosaminyl-transferase V antibody was raised against a peptide at the C terminal residues of human N-acetylglucosaminyl-transferase V protein.
Target Name	N-acetylglucosaminyltransferase V
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

## Application Details

ELISA IHC

Note: This product is for in vitro research use only