

hsp60 Antibody

Catalog No: #39380

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	hsp60 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	WB IHC
Species Reactivity	Hu
Specificity	Human
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide
Target Name	hsp60
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded breast carcinoma showing cytoplasmic staining with anti-Hsp60 antibody at a dilution of 1/100(A). Block with immunizing peptide(B).



Note: This product is for in vitro research use only