

# SCP2 Antibody

Catalog No: #39403

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	SCP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	WB IHC
Species Reactivity	Hu
Specificity	Human
Target Name	SCP2
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

## Application Details

WB IHC

## Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal testis showing nuclear with anti-SCP2 antibody at a dilution of 1/100.

### HepG2 2.2.15 Fetal testis



Western Blot analysis of HepG2 2.2.15 cell and fetal testis Lysate with anti-SCP2 antibody

WB 1:500

---

Note: This product is for in vitro research use only