

## RBBP6 Antibody

Catalog No: #39408

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	RBBP6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	IHC
Species Reactivity	Hu Ms
Specificity	Human, Mouse
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein
Target Name	RBBP6
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

## Application Details

IHC

## Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal testis showing nuclear staining with anti-RBBP6 antibody at a dilution of 1/100.

Note: This product is for in vitro research use only