

## ABCD3 Antibody

Catalog No: #39411

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ABCD3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	WB IHC
Species Reactivity	Hu
Specificity	Human
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein
Target Name	ABCD3
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


## Application Details

WB IHC

## Images



**Immunohistochemical staining of formalin-fixed paraffin-embedded fetal liver showing membrane staining with anti-ABCD3 antibody at a dilution of 1/100.**



---

Note: This product is for in vitro research use only