

myocardin Antibody

Catalog No: #39414

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	myocardin Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	MYCD Polyclonal Antibody detects endogenous levels of protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide
Target Name	myocardin
Uniprot	Q8IZQ8
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at +4°C for short term. Store at -20°C for long term. Avoid freeze/thaw cycle.

Application Details

WB 1:500-2000

Images



Western Blot analysis of Hela lysis, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only