

## 2-hydroxyacyl-CoA lyase 1 Antibody

Catalog No: #39416

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	2-hydroxyacyl-CoA lyase 1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	WB IHC
Species Reactivity	Hu Ms
Specificity	Human?????B'B—Mouse
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein
Target Name	2-hydroxyacyl-CoA lyase 1
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


## Application Details

WB IHC

## Images



**Immunohistochemical staining of formalin-fixed paraffin-embedded fetal liver showing cytoplasmic staining with anti-2 Hydroxy phytanoyl CoA lyase antibody at a dilution of 1/100.**



---

Note: This product is for in vitro research use only