

GALNT9 Antibody

Catalog No: #39419

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GALNT9 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	WB IHC
Species Reactivity	Hu Ms
Specificity	Human
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein
Target Name	GALNT9
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details


WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal brain showing Cytoplasmic and Membrane staining with anti-GALNT9 antibody at a dilution of 1/100.

Mouse Fetal
brain brain



— 95KD
— 72KD
— 55KD
— 43KD

Western Blot analysis
of mouse brain and
fetal brain Lysate with
anti-GALNT9 antibody

WB 1:1000

Note: This product is for in vitro research use only