

phospholipid scramblase 1 antibody

Catalog No: #22435

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	phospholipid scramblase 1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 256 of phospholipid scramblase 1
Target Name	phospholipid scramblase 1
Accession No.	Swiss-Prot:O15162Gene ID:5359
Uniprot	O15162
GeneID	5359;
Concentration	0.5mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details


Predicted MW: 35kd

Western blotting: 1:500-1:3000


Immunohistochemistry: 1:100-1:500

Immunofluorescence: 1:100-1:200


Images



Sample(30 ug whole cell lysate)
A: HeLa S3
10% SDS PAGE
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded MDAMB468 xenograft, using phospholipid scramblase 1 antibody at 1: 500 dilution.



Immunofluorescence analysis of methanol-fixed A431, using Scramblase1 antibody at 1: 200 dilution.

Note: This product is for in vitro research use only