

# HHLA3 Antibody

Catalog No: #39441

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	HHLA3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu
Specificity	Human
Target Name	HHLA3
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

## Application Details

ELISA WB IHC

## Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal kidney showing cytoplasmic and nuclear staining with anti-HHLA3 antibody at a dilution of 1/100.

fetal fetal fetal  
kidney liver heart



Western Blot  
analysis of fetal  
kidney, liver and  
heart Lysate with  
anti-HHLA3  
antibody

WB 1:1000

---

Note: This product is for in vitro research use only