

## HERC6 Antibody

Catalog No: #39442

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	HERC6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified (Protein A purification)
Applications	ELISA WB IHC
Species Reactivity	Hu
Specificity	Human
Target Name	HERC6
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


## Application Details

ELISA WB IHC

## Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal testis showing cytoplasmic staining with anti-HERC6 antibody at a dilution of 1/100.



---

Note: This product is for in vitro research use only