

ADORA2A Antibody

Catalog No: #39487

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ADORA2A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu
Specificity	ELISA titerWestern blot analysisImmunohistochemistry
Target Name	ADORA2A
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA: 1:20000-1:80000 WB: 1:500-1:1000 IHC: 1:100-1:500

Note: This product is for in vitro research use only