


## Description

Product Name	AGAP3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Rt Ms
Target Name	AGAP3
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

## Application Details


ELISA WB IHC

## Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal brain showing cytoplasmic staining with anti-AGAP3 antibody at a dilution of 1/100.

fetal fetal fetal  
heart lung kidney



Western Blot analysis of fetal heart, lung and kidney Lysate with anti-AGAP3 antibody

WB 1:500

---

Note: This product is for in vitro research use only