

PCBP1 Antibody

Catalog No: #39508

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PCBP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Ms Rt
Target Name	PCBP1
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal thymus showing cytoplasmic and nucleus staining with anti-PCBP1 antibody at a dilution of 1/100.

Hela K-562



WB 1:500

Western Blot analysis of Hela and K-562 cell Lysate with anti-PCBP1 antibody

Note: This product is for in vitro research use only