

SCFD1 Antibody

Catalog No: #39510

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SCFD1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Ms Rt
Target Name	SCFD1
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details


ELISA WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal pancreas showing cytoplasmic staining with anti-SCFD1 antibody at a dilution of 1/100.

Hela HepG2 293T



Western Blot analysis of Hela, HepG2 and 293T cell Lysate with anti-SCFD1 antibody

WB 1:500

Note: This product is for in vitro research use only