

POLR3B Antibody

Catalog No: #39526

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	POLR3B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Ms Rt
Target Name	POLR3B
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details


ELISA WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal liver showing nucleus staining with anti-POLR3B antibody at a dilution of 1/100.

293T fetal liver



Western Blot analysis of 293T cell and fetal liver Lysate with anti-POLR3B antibody

WB 1:500

Note: This product is for in vitro research use only