

WNK1 Antibody

Catalog No: #39539

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	WNK1 Antibody
Host Species	Rabbit
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant human WNK1 protein 80-338aa (BC130467)
Target Name	WNK1
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details


ELISA WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal kidney showing cytoplasmic staining with anti-WNK1 antibody at a dilution of 1/100.

HEK293 293T HepG2



Western Blot analysis of HEK293, 293T and HepG2 cell Lysate with anti-WNK1 antibody

WB 1:500

Note: This product is for in vitro research use only