

PELI3 Antibody

Catalog No: #39541

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PELI3 Antibody
Host Species	Rabbit
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant human PELI3 protein 52-299aa (BC021256)
Target Name	PELI3
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details


ELISA WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal heart staining with anti-PELI3 antibody at a dilution of 1/100.

MCF-7 fetal fetal
thyroid colon



Western Blot
analysis of MCF-7
cell, fetal thyroid
and colon Lysate
with anti-PELI3
antibody

WB 1:500

Note: This product is for in vitro research use only