

IMP4 Antibody

Catalog No: #39551

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IMP4 Antibody
Host Species	Rabbit
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant Protein
Immunogen Description	Raised against a recombinant human IMP4 protein 102-281aa (BC010042)
Target Name	IMP4
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

ELISA WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal spleen showing nucleus staining with anti-IMP4 antibody at a dilution of 1/100.



Note: This product is for in vitro research use only