

INPP4A Antibody

Catalog No: #39554

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	INPP4A Antibody
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant Protein
Immunogen Description	Raised against a recombinant human INPP4A protein 601-812aa (BC028361)
Target Name	INPP4A
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

ELISA WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal lymphocyte staining with anti-INPP4A antibody at a dilution of 1/100.



Note: This product is for in vitro research use only