

## PDK3 Antibody

Catalog No: #39575

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	PDK3 Antibody
Host Species	Rabbit
Purification	Affinity purified
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant Protein
Immunogen Description	Raised against a recombinant human PDK3 protein 102-351aa (BC015948).
Target Name	PDK3
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

## Application Details

WB IHC

## Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal skeletal muscle showing cytoplasmic staining with anti-PDK3 antibody at a dilution of 1/100.

HepG2



-95KD  
-72KD  
-55KD  
-43KD

Western Blot analysis of  
HepG2 cell Lysate with  
anti-PDK3 antibody

WB 1:1000

---

Note: This product is for in vitro research use only