

IL34 Antibody

Catalog No: #39576

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IL34 Antibody
Host Species	Rabbit
Purification	Affinity purified
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant Protein
Immunogen Description	Raised against a recombinant human IL34 protein 21-242aa (BC029804).
Target Name	IL34
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal spleen staining with anti-IL34 antibody at a dilution of 1/100.

fetal brain



Western Blot analysis of fetal brain Lysate with anti-IL34 antibody

WB 1:500

Note: This product is for in vitro research use only