

## IFIT2 Antibody

Catalog No: #39584

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	IFIT2 Antibody
Purification	Affinity purified
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant Protein
Immunogen Description	Raised against a recombinant human IFIT2 protein 100-321aa (BC032839).
Target Name	IFIT2
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

## Application Details

WB: 1:200~1:1000

IHC: 1:100-1:500

## Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal liver showing cytoplasmic staining with anti-IFIT2 antibody at a dilution of 1/100.



Western Blot analysis of fetal liver and lung Lysate with anti-IFIT2 antibody

WB 1:1000

---

Note: This product is for in vitro research use only