

# ITIH2 Antibody

Catalog No: #39623

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ITIH2 Antibody
Host Species	Rabbit
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant Protein
Immunogen Description	Raised against a recombinant human ITIH2 protein 124-321aa (BC132685).
Target Name	ITIH2
Concentration	1.0mg/ml
Formulation	Reconstitute in 100ul of sterile H <sub>2</sub> O with 50% glycerol 0.01% thiomersal.
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles


## Application Details

WB 1:500-1:2,000


IHC 1:100-1:200

ELISA 1:20,000-1:80,000

## Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal brain staining with anti-ITIH2 antibody at a dilution of 1/100.



---

Note: This product is for in vitro research use only