

## PIH1D2 Antibody

Catalog No: #39639

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	PIH1D2 Antibody
Host Species	Rabbit
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant Protein
Immunogen Description	Raised against a recombinant human PIH1D2 protein 25-313aa (BC019238).
Target Name	PIH1D2
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


## Application Details

ELISA WB IHC

## Images



**Immunohistochemical staining of formalin-fixed paraffin-embedded fetal testis staining with anti-PIH1D2 antibody at a dilution of 1/100.**



---

Note: This product is for in vitro research use only