

ZFP28 Antibody

Catalog No: #39694

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZFP28 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Rt Ms
Immunogen Type	Recombinant Protein
Immunogen Description	raised against a recombinant human ZFP28 protein 21-278aa (BC028416).
Target Name	ZFP28
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

ELISA WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal colon showing nuclear staining with anti-ZFP28 antibody at a dilution of 1/100.



WB: 1:1,000

Note: This product is for in vitro research use only