

SZT2 Anibody

Catalog No: #39697

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SZT2 Anibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Rt Ms
Immunogen Type	Recombinant Protein
Immunogen Description	raised against a recombinant human SZT2 protein 299-479aa (BC052802).
Target Name	SZT2 Anibody
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

ELISA WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal brain showing cytoplasmic staining with anti-SZT2 antibody at a dilution of 1/100.



Note: This product is for in vitro research use only