

## KIAA0284 Antibody

Catalog No: #39700

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	KIAA0284 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant Protein
Immunogen Description	raised against a recombinant human KIAA0284 protein 1444-1583aa (BC047913)
Target Name	KIAA0284
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


## Application Details

ELISA WB IHC

## Images



**Immunohistochemical staining of formalin-fixed paraffin-embedded fetal skeletal muscle showing cytoplasmic staining with anti-KIAA0284 antibody at a dilution of 1/100.**



---

Note: This product is for in vitro research use only