

KHDC1 Antibody

Catalog No: #39701

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	KHDC1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant Protein
Immunogen Description	raised against a recombinant human KHDC1 protein 28-158aa (BC022080).
Target Name	KHDC1
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

ELISA WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal kidney staining with anti-KHDC1 antibody at a dilution of 1/100.



Note: This product is for in vitro research use only