

PISD Antibody

Catalog No: #39713

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PISD Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen affinity column
Applications	ELISA WB IHC
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant Protein
Immunogen Description	raised against a recombinant human PISD protein 71-310aa (BC001482)
Target Name	PISD
Concentration	1ug/ul
Formulation	Reconstitute in 100ul of sterile H ₂ O with 50% glycerol 0.01% thiomersal.
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

WB α O 1:500-1:2,000IHC α O 1:100-1:200

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal kidney showing cytoplasmic staining with anti-PISD antibody at a dilution of 1/100.



Note: This product is for in vitro research use only