

OX40 Antibody

Catalog No: #39734

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description


Product Name	OX40 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB
Specificity	Human and predicted against mouse and rat.
Immunogen Type	Peptide-KLH
Immunogen Description	Rabbit polyclonal OX40 (1) antibody was raised against a peptide from C terminal residues of human OX40 protein.
Target Name	OX40
Uniprot	P43489
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA: 1:20000-1:80000 WB: 1:500-1:1000 IHC^o 1:100-1:200

Images





Immunohistochemical staining of formalin-fixed paraffin-embedded fetal spleen showing membrane staining

Note: This product is for in vitro research use only