

AMN1 Antibody

Catalog No: #39755

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	AMN1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA IHC
Specificity	Human and predicted against mouse and rat.
Immunogen Type	Peptide-KLH
Immunogen Description	Rabbit polyclonal AMN1 (1) antibody was raised against a peptide near the intermediate corresponding to the C terminus residues of human AMN1.
Target Name	AMN1
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA: 1:20000-1:80000 IHC: 1:100-1:500

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded colon carcinoma showing membrane staining with anti-AMN1 antibody at a dilution of 1/100.

Note: This product is for in vitro research use only