

KIAA0802 Antibody

Catalog No: #39785


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	KIAA0802 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Specificity	Human and predicted against mouse and rat.
Immunogen Type	Recombinant Protein
Immunogen Description	Rabbit polyclonal KIAA0802 (1) antibody was raised against a recombinant human KIAA0802 protein 13-202aa (BC040542).
Target Name	KIAA0802
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal testis staining with anti-KIAA0802 antibody at a dilution of 1/100



Note: This product is for in vitro research use only