

## USP29 Antibody

Catalog No: #39787

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	USP29 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Specificity	Human and predicted against mouse and rat.
Immunogen Type	Recombinant Protein
Immunogen Description	Rabbit polyclonal USP29 (1) antibody was raised against a recombinant human USP29 protein 15-215aa (BC130394).
Target Name	USP29
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


## Application Details

ELISA titer: 1:100


Western blot analysis: 1:1,000

Immunohistochemistry: 1:100 optimal dilutions/concentrations should be determined by the end user.

## Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal testis staining with anti-USP29 antibody at a dilution of 1/100



---

Note: This product is for in vitro research use only