

LENG8 Antibody

Catalog No: #39875

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LENG8 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen affinity column
Applications	ELISA WB IHC
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant Protein
Immunogen Description	Rabbit polyclonal LENG8 (1) antibody was raised against a recombinated human LENG8 protein 489-664aa (BC028048).
Target Name	LENG8
Concentration	1.0mg/ml
Formulation	Reconstitute in 100ul of sterile H2O with 50% glyceroloO 0.01% thiomersal.
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles


Application Details

WB 1:500-1:2000



IHC 1:100-1:200

ELISA 1:20,000-1:80000

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal kidney showing cytoplasmic staining with anti-LENG8 antibody at a dilution of 1/100



Note: This product is for in vitro research use only