

ZNF276 Antibody

Catalog No: #39884

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZNF276 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen affinity column
Applications	ELISA WB IHC
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant Protein
Immunogen Description	Rabbit polyclonal ZNF276 (1) antibody was raised against a recombinated human ZNF276 protein 175-385aa (BC107603).
Target Name	ZNF276
Calculated MW	60kDa
Concentration	1.0mg/ml
Formulation	Reconstitute in 100uL of sterile H ₂ O with 50% glycerol/0.01% thimerosal.
Storage	For extended storage, aliquot contents and freeze at -20/ C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

ELISA: 1:20000-1:80000;


WB: 1:500-1:2000;

IHC: 1:100-1:200

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal tonsil showing nuclear staining with anti-ZNF276 antibody at a dilution of 1/100



Note: This product is for in vitro research use only