

## PPP1R3C Antibody

Catalog No: #39988

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	PPP1R3C Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen affinity column
Applications	ELISA WB IHC
Species Reactivity	Hu Rt Ms
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fragment corresponding to 79-295aa of human PPP1R3C protein
Target Name	PPP1R3C
Concentration	1.0mg/ml
Formulation	Reconstitute in 100ul of sterile H <sub>2</sub> O with 50% glycerol 0.01% thiomersal.
Storage	Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles



## Application Details

WB 1:500-1:2000

IHC 1:100-1:200

ELISA 1:20000-1:80000

## Images



---

Note: This product is for in vitro research use only